

Product Specifications



Specification	Detail
Product Code	GOB08
Product Name	Surgical Laparotomy Sponges
Material	100% Prewashed Cotton
Size	18" x 18"
Sterility	Sterile, Single-Use
X-Ray Detectable	Yes, Blue Radiopaque Thread
Loop	7-inch Cotton Loop
Color	White with Blue Indicator
Packaging	5 pcs / Sterile Soft Pack
Inner Box	40 Bags / Blister Printed Bag
Master Carton	200 piece (40 bags) / Carton
Latex-Free	Yes
Shelf Life	5 Years
Manufactured For	Goben Medical, Aventura, FL 33160, USA

R_x ONLY

Prescription
use only



Single-use
only



Do not
reprocess



Do not use if
packaging is
damaged



Sterile, ethylene
oxide

