

Product Specifications



Specification	Detail
SKU	GOB03
Type	Medical Suction Tubing
Configuration	Male-to-Male Connectors
Inner Diameter	1/4 inch (6.35 mm)
Length	10 ft (3.05 m)
Material	Medical-grade Clear PVC
Sterilization	EO (Ethylene Oxide) Sterile
Use	Single-use / Disposable
Packaging	1 pc per sterile blister, 50 pcs per carton
Color	Clear Transparent
Compatibility	Yankauer Tips, Suction Canisters, Vacuum Regulators
Compliance	Manufactured under ISO / FDA standards
Shelf Life	5 years
Manufactured For	Goben Medical, Aventura, FL 33160, USA



Prescription
use only



Non-pyrogenic



Single-use
only



Do not
reprocess



Do not use if
packaging is
damaged



Sterile, ethylene
oxide



GOBEN
MEDICAL