

Product Specifications



Specification	Detail
Part Number	GOB01
Type	Infiltration Tubing for Peristaltic Pump
Configuration	Single Spike, Luer Lock Fitting
Length	9 ft (2.74 m)
Material	Medical-grade, kink-resistant tubing
Sterilization	EO (Ethylene Oxide) Sterile
Use	Single Use / Disposable
Packaging	1 pc per sterile blister pack
Inner Box	10 pcs per white box
Master Carton	100 pcs per carton
Manufactured For	Goben Medical, Aventura, FL 33160, USA

R_x ONLY

Prescription
use only



Non-pyrogenic



Single-use
only



Do not
reprocess



Do not use if
packaging is
damaged



Sterile, ethylene
oxide



GOBEN
MEDICAL