

Product Specifications



Specification	Detail
Product Code	GOB09
Description	Surgical Gown AAMI Level 3, Non-Reinforced, with Towel
Material	50 g SMS Non-Woven Polypropylene
Size	52 in – Extra Large (XL)
Color	Blue
Protection Level	AAMI Level 3 (ANSI/AAMI PB70)
Design	Non-Reinforced, Set-in Sleeves, Knit Cuffs
Closure	Adjustable Neck Hook-and-Loop + Waist Ties
Included Accessories	1 Towel (per pack)
Sterility	Sterile – Ethylene Oxide (EO)
Packaging	1 pc / Sterile Blister Pack with Label
Carton Quantity	40 pcs / Carton
Latex-Free	Yes
Compliance	ISO 13485, CE, FDA Registered Manufacturing
Shelf Life	5 Years
Manufactured For	Goben Medical, Aventura, FL 33160, USA



Prescription
use only



Non-pyrogenic



Single-use
only



Do not
reprocess



Do not use if
packaging is
damaged



Sterile, ethylene
oxide



GOBEN
MEDICAL