

Product Specifications



Specification	Detail
Product Code	GOB05
Product Type	Endotracheal Tube (Cuffed)
Inner Diameter (I.D.)	6.5 mm
Material	Medical-grade PVC, latex-free
Cuff Type	High volume / low pressure
Tip Type	Murphy eye
Curve Type	Magill curve
Pilot Balloon	Advanced self-sealing valve
Radiopaque Line	Yes
Markings	Centimeter depth graduations
Use	Oral or nasal intubation
Sterilization Method	Ethylene Oxide (EO)
Usage	Single-use, sterile, disposable
Packaging	1 pc / printed blister bag
Inner Box	10 pcs / white box with label
Master Carton	100 pcs / carton
Shelf Life	5 years
Manufactured For	Goben Medical, Aventura, FL 33160, USA

Rx ONLY

Prescription
use only



Single-use
only



Do not
reprocess



Do not use if
packaging is
damaged



Sterile, ethylene
oxide

