

# Product Specifications



Specification	Detail
<b>Part Number</b>	GOB04
<b>Product Type</b>	Closed Suction Wound Drainage Kit
<b>Reservoir Capacity</b>	100ml
<b>Drain Type</b>	Flat 10mm, 3/4 perforations
<b>Material</b>	Medical-grade silicone, latex-free
<b>Sterilization</b>	Ethylene Oxide (EO)
<b>Usage</b>	Single-use, sterile, disposable
<b>Packaging</b>	1 pc / printed blister pack
<b>Inner Box</b>	10 pcs / white box with label
<b>Master Carton</b>	100 pcs / carton
<b>Comparison</b>	Comparable to Jackson-Pratt® style systems
<b>Compliance</b>	ISO 13485 / FDA-registered manufacturing
<b>Shelf Life</b>	5 years
<b>Manufactured For</b>	Goben Medical, Aventura, FL 33160, USA

**Rx ONLY**

Prescription use only



Single-use only



Do not reprocess



Do not use if packaging is damaged



Sterile, ethylene oxide