

Product Specifications



Specification	Detail
Part Number	GOB02
Length	9 ft
Internal Diameter	3/8 in
Material	Super soft, medical-grade, kink-resistant tubing
Sterility	Sterile, single-use
Latex	Latex-Free
Memory	Reduced "memory" for smooth handling
Packaging	1 pc per printed blister pack
Master Carton	20 pcs per carton
Compatibility	Compatible with all standard liposuction aspirators
Recommended Use	Liposuction, tumescent infiltration, and other high-vacuum procedures
Shelf Life	5 years
Manufactured For	Goben Medical, Aventura, FL 33160, USA

R_x ONLY

Prescription
use only



Non-pyrogenic



Single-use
only



Do not
reprocess



Do not use if
packaging is
damaged



Sterile, ethylene
oxide



GOBEN
MEDICAL