

# Product Specifications



Specification	Detail
<b>Product Code</b>	GOB09
<b>Description</b>	Surgical Gown AAMI Level 3, Non-Reinforced, with Towel
<b>Material</b>	50 g SMS Non-Woven Polypropylene
<b>Size</b>	52 in – Extra Large (XL)
<b>Color</b>	Blue
<b>Protection Level</b>	AAMI Level 3 (ANSI/AAMI PB70)
<b>Design</b>	Non-Reinforced, Set-in Sleeves, Knit Cuffs
<b>Closure</b>	Adjustable Neck Hook-and-Loop + Waist Ties
<b>Included Accessories</b>	1 Towel (per pack)
<b>Sterility</b>	Sterile – Ethylene Oxide (EO)
<b>Packaging</b>	1 pc / Sterile Blister Pack with Label
<b>Carton Quantity</b>	50 pcs / Carton
<b>Latex-Free</b>	Yes
<b>Compliance</b>	ISO 13485, CE, FDA Registered Manufacturing
<b>Shelf Life</b>	5 Years
<b>Manufactured For</b>	Goben Medical, Aventura, FL 33160, USA



Prescription  
use only



Non-pyrogenic



Single-use  
only



Do not  
reprocess



Do not use if  
packaging is  
damaged



Sterile, ethylene  
oxide



**GOBEN**  
MEDICAL