

Product Specifications



Specification	Detail
Product Code	GOB06
Product Name	Guedel Airway, Oropharyngeal
Size	90mm (Adult, ISO Size 9)
Color Code	Yellow
Material	Medical-Grade Plastic (Latex-Free)
Type	Disposable, Non-Sterile
Functions	Airway, Bite Block, Tongue Depressor
Center Channel	Yes, for suction catheter passage
Packaging	1 pc / PE pack
Inner Box	50 pcs / White box with label
Master Carton	250 pcs / Carton
Shelf Life	5 Years
Compliance	ISO 13485, CE Certified, FDA Registered Manufacturing
Manufactured For	Goben Medical, Aventura, FL 33160, USA

Rx ONLY

Prescription
use only



Single-use
only



Do not
reprocess



Do not use if
packaging is
damaged



Sterile, ethylene
oxide



GOBEN
MEDICAL