


Product Specifications



Specification	Detail
Product Code	GOB10
Material	PE Film 50g + SMS Non-Woven 33g
Dimensions	44 in x 90 in
Color	Blue
Sterility	EO Sterile
Reinforcement Type	Zone Reinforced
Usage	Disposable, Single Use
Packaging	1 pc / Blister Pack
Carton Quantity	25 pcs / ctn
Latex Free	Yes
Shelf Life	5 Years
Manufactured For	Goben Medical, Aventura, FL 33160, USA

R_x ONLY

Prescription
use only



Non-pyrogenic



Single-use
only



Do not
reprocess



Do not use if
packaging is
damaged



Sterile, ethylene
oxide



GOBEN
MEDICAL